




InfantRisk Center | HCP App Screenshots


Home Screen



This is the home screen of the application. From here you can navigate to a keyword search for prescription and non-prescription products, patient customization, a direct dial to the InfantRisk Center, a list of hot topics, your favorite bookmarks, and your profile.



Search



From the search screen, you can either search by keyword or medication name, or browse all listings on prescription drugs, non-prescription drugs, or vitamins. Additionally, you can search by condition to see listings that fall under specific ailments.



Search Results by Patient

AT&T 3G 1:48 PM

Search Prescription Drugs ⓘ

You are looking at information that is customized for Pam | [New Patient Search](#)

- D-50
- DABIGATRAN ETEXILATE
- Dacodyl
- DACTINOMYCIN
- DALFOPRISTIN
- Dalmane

Home Search Patient Call IRC More

This is the search results for a specific patient's situation. Results are color coded with safety ratings on the left, based on the safety of the product for the specific patient.



AT&T 3G 10:32 AM

Results Details

1ST 2ND 3RD BIRTH 0-6 6-12 12 +
TRIMESTER MONTHS


DRUG NAME: IBUPROFEN
DRUG TYPE: Analgesic, antipyretic
TRADE
NAMES: Advil, Nuprin, Motrin, Pediafen
USUAL DOSE: 400 mg every 4-6 hours
LACTATION RISK: L1
RELATIVE INFANT DOSE: 0.1181% - 0.6562%
SIDE EFFECTS: Nausea, epigastric pain, dizziness, edema, GI bleeding.
ALTERNATIVE MEDS: Acetaminophen
T1/2: 1.8-2.5 hours
ORAL BIOAVAILABILITY: 80%
MW: 206

Home Search Patient Call IRC More

This is the drug detail screen. It shows safety ratings for each stage of pregnancy and nursing, along with descriptions of the drug, dosage information, side effects and alternative options.



Hot Topics



The Hot Topics page displays the latest Hot Topics articles posted to InfantRisk.com, and displays a mobile version of the articles on the site when clicked.



Call the IRC



You can also call the InfantRisk center. Allows the user to instantly call the IRC at the click of the button, and displays the local time at the center, so users can easily determine whether or not it is open.